LUMATEC

---

SUPERLITE  
M 05

---

THE SMALLEST  
MULTISPECTRAL LED  
LIGHT SOURCE

Operating Instructions

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[EN] These operating instructions belong with this product. They contain important information for putting it into service and operating it. This should be noted also when this product is passed on to a third party.

Therefore look after these operating instructions for future reference!

A list of contents with the corresponding page numbers can be found in the index on page 1.

LUMATEC prepared these instructions for experienced forensic scientists, process engineers, technicians and manufacturing personnel. If you are new to UV light applications, do not understand the instructions or have any further questions and concerns, you can find our contact information on the last page of this booklet.

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**Part No**

2350.0017	SUPERLITE M 05	[Power Pack I]
2350.0016	SUPERLITE M 05	[Charge Adapter and Charger for Power Pack I]
2350.0060	SUPERLITE M 05	[Power Pack II]
2350.0061	SUPERLITE M 05	[Charger for Power Pack II]
2350.00XX	SUPERLITE M 05	[LED Head]

---

## TABLE OF CONTENTS

1.	Intended Use.....	2
2.	Safety Instructions.....	2
3.	Equipment Overview .....	3
4.	Equipment Installation .....	5
5.	Equipment Operation .....	7
6.	Accessories and Spare Parts .....	9
7.	Information on UV Radiation.....	10
8.	Spectral Ranges and Applications.....	11
9.	Technical Data.....	12

---

## 1. INTENDED USE

The light source SUPERLITE M 05 is solely intended for fluorescence excitation in laboratories or at crime scenes and for technical inspection as well as curing of polymers in industrial environments.

---

## 2. SAFETY INSTRUCTIONS

Before using the SUPERLITE M 05, please read the following chapters of safety, installation and operation instructions.

### Safety Symbols

Warning!

Potential Danger

Warning!

Optical Radiation

Mandatory Action!

Wear UV protection goggles

### Mains Supply

The intelligent electronic power supply of the unit will adapt automatically to all AC currents from 110 to 240 Volts and 50 to 60 Hz.

### Explosive Surroundings

This unit is not meant for operation in explosive surroundings.

### Humid Surroundings

This unit is only intended for use in dry environments.

### Protection against Strong Optical Radiation

Never look directly into the light. Shield the operator or the light source. To avoid tissue damage, do not expose the unprotected eye or skin to the ultraviolet light. When work under unshielded UV radiation is necessary, wear UV protection goggles, long-sleeved shirt or lab coat and gloves.

When using the 405 nm and 455 nm lens heads (violet and blue), we strongly recommend wearing yellow or orange filter goggles. By suppressing the excitation wavelength, fluorescence is perceived much better due to increased contrast and at the same time the user is protected against unintentional glare.

### Danger of Fire

Do not deposit on inflammable objects like cloth or paper. The intensive radiation may be sufficient to set these objects on fire.

### 3. EQUIPMENT OVERVIEW


#### Nomenclature

1	Lense Ring	6	Ventilation Exhaust
2	Wavelength Selector	7	1/4" UNC threaded Hole
3	Anti-Twist Pin	8	Status LED
4	Loop for Lanyard	9	Charging / Mains Socket
5	Push-Button	10	New hexagonal Design

#### Front and Back View

##### LED-Single-Wavelength-Head

For only one wavelength with greater intensity


##### LED-Multihead


For 4, 5 or 6 different wavelengths

##### Power Pack I

##### Power Pack II

Whether you are looking for forensic evidence, performing leak tests or conducting nondestructive tests – the SUPERLITE M 05 is the only one in its performance class to feature an interchangeable head system that allows you to quickly and easily choose between different spectra. This state-of-the-art LED light source is equipped with high-quality optics to emit homogeneous light. The SUPERLITE M 05 offers three intensity levels, so you can adjust the light's intensity to suit the given requirements.

### LED-Heads

Available are LED-Multiheads that contain LEDs for 4, 5 or 6 different wavelengths and LED-Single-Wavelength-Heads that produce radiation of only one wavelength, but with greater intensity.

### Power Pack - Version I and II

The Power Pack also contains the control electronics for the LED-Head, a cooling fan and the Push-Button to switch the lamp on and off and to regulate the intensity in three levels. While the old version I is charged using the Charger Adapter, this is already integrated in the new version.

Other advantages of version II are:

- New hexagonal design (10) to prevent rolling on flat surfaces.
- Suitable for continuous operation with the help of the charging / mains connection socket (9).
- Charging status can be checked via LED (8).
- Integrated 1/4" UNC threaded hole (7) for mounting on a tripod
- The Anti-Twist Pin (3) enables a twist-proof connection to the LED-Multiheads, which effectively prevents damage due to incorrect handling.

Both Power Packs (versions I and II) can be combined with a LED-Single-Wavelength-Head or a LED-Multihead.

**Note:** If you want to combine the newer Power Pack II with an older Single-Wavelength-Head, the Anti-Twist Pin (3) may have to be unscrewed from the housing of the Power Pack II, as older Single-Wavelength-Heads did not have the recess for the pin.

## 4. EQUIPMENT INSTALLATION

### Charging the Power Pack I

On delivery the Power Pack I is only partially charged. Before first use it must be charged with the Charger Adapter and the AC/DC Power Supply.


The three contact pins of the Charger Adapter fit into the three sockets of the Power Pack I and snap in with a clicking sound. Align the two engraved arrows on both parts to find the correct position.


Plug the Charger Adapter into the Power Pack I and connect it to the round plug of the AC/DC Power Supply (charge voltage 5V!).

With the SUPERLITE M regulated Car Charger Adapter it can be charged in a car with engine off.

The flashing indicator LED of the Charger Adapter shows the Power Pack I is being charged. The Power Pack I is fully charged when the LED is continuously on.


### Connecting the LED-Heads

The LED-Heads are connected to the Power Pack I in the same way as the Charger Adapter.

**Important:** To connect or disconnect, always push or pull, never twist!

After connecting the Power Pack I to a LED-Head, the SUPERLITE M 05 is ready to be switched on.

## Charging the **Power Pack II**

On delivery the Power Pack II is only partially charged. Before first use it must be charged with the help of the AC/DC Power Supply, which is also called "Charger". Plug the round connector of the AC/DC Power Supply (charge voltage 5V!) into the Charging Socket (9) of the Power Pack II. The flashing Status LED (8) shows, that the Power Pack II is being charged. The Power Pack II is fully charged when the LED is continuously on.

An LED-Lens-Head may or may not be installed during charging. While charging, the fan is running to cool the unit.

## Mains operation with **Power Pack II**

If an LED-Lens-Head is installed during charging and the lamp is switched on by using the Push-Button (5), the Power Pack II switches to mains operation.

The charging process is delayed. The fan runs at different speeds depending on the temperature.


## Installing an LED-Lens-Head (Single-Wavelength or Multi-head) on the **Power Pack II**

Align the two arrows engraved on the housing and slide the LED-Lens-Head onto the Power Pack II. In this way, the three pins of the head are guided into the three sockets of the Power Pack, thus establishing the electrical connection. In case connecting is not possible, please observe the note on page 4 about the Anti-Twist Pin (3).

To disconnect, pull both units apart in the direction of the arrows.

**Important:** To connect or disconnect, always push or pull, never twist!

**Note:** The Charger Adapter of Power Pack I has no function when used with Power Pack II.

## 5. EQUIPMENT OPERATION

The SUPERLITE M 05 can be switched on by pressing the Push-Button.

**Note:** Before switching on the LED-Multihead, make sure that a locking position is reached. If this is not the case, the electrical circuit is not closed and the lamp cannot be switched on.

The Push-Button reacts to

- A short tap
- Pressing and holding the Push-Button

Activating the LED-Heads

Light is off


- 1 x tap Push-Button: 25% of luminous flux, fan turns on

Light is on (25% or 50% Mode)

- 1 x tap Push-Button: Light turns off, fan turns off
- 1 x press and hold Push-Button: Light switches into the next higher power mode (50%, 100%)

Light is in 100% Mode

- Without action: After 60 seconds light switches to 25% mode
- 1 x tap Push-Button: Light turns off, fan turns off
- 1 x press and hold Push-Button: Light switches to 25% mode

## LED-Single-Wavelength-Heads, Optical Power Output, Runtime

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
## LED-Multiheads for several Wavelengths, Optical Power Output, Runtime

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### Drop and Roll Safety

The lanyard which can be hooked to the Power Pack I with the help of its loop (4) is meant to be worn around the neck or the wrist. Never remove this lanyard. It prevents the light source from rolling off inclined surfaces. The Power Pack II has a hexagonal design which prevents rolling, here the lanyard is not necessary.

### Outdoor Use

If water entered the air vents during use, the light source should be left on to dry inner components. Finally the LED-Head and Power Pack should be detached, shaken out by hand and dried in warm air to avoid corrosion.

### Incident and Line-Illumination for Forensics

The coloured Lense Ring (1) of a LED-Single-Wavelength-Head can be pushed back and forward by 10 mm to focus the light. This way the light beam can be adjusted from a narrow angle (pushed in) to a wide homogeneous beam (pulled out).

**Note:** LED-Multiheads do not have this feature.

Even tiny objects like hair, skin particles, ridges of fingerprints in blood, fat or other consolidating substances, shoe or tire prints can be visualized and photographed.

For detailed inspections or photography, use the white LED-Head (400-700 nm) with up to 100% power and pulled out Lense Ring for homogeneous illumination.

### Scanning of Larger Areas

When searching objects at a long distances, the white LED-Head (400-700 nm) with pushed-in Lense Ring may be used in 100% Mode, in order to work with a far reaching beam which is as powerful as possible.

An optional “Ultrathrow” LED-Head is available, optimized for long distances (search light).

### Setting the Wavelength on LED-Multiheads

Depending on the model, the LED-Multihead includes 4, 5 or 6 different LEDs with different wavelengths, integrated in one head. By simply turning the housing, the excitation wavelength can be changed.

**Note:** If you turn the housing very slowly between two stops, the head will automatically switch off for safety reasons to protect electrical components. By turning the housing to a mark and pressing the Push-Button on the Power Pack, the lamp can be switched on again.

**Note:** Before switching on the LED-Multihead, make sure that a locking position is reached. If this is not the case, the electrical circuit is not closed and the lamp cannot be switched on.

### Fluorescence and Contrast Generation

Fluorescence signals are generated with wavelengths of 365 nm, 390 nm (UVA), 405 nm (violet) and 445 nm (blue).

Biological fluorescence signals are generated by illumination with 445 nm (blue) and observed with a pair of orange filter goggles. The evidence is logged by camera through an orange camera filter. The intensity of the excitation is raised by switching to 50% or 100% Mode and additionally pushing in the Lense Ring if applicable.

Many objects show good contrast against fluorescing backgrounds or are visualized as black spots by illumination with

a colour according to their chemical absorption maximum (e.g. blood stains illuminated with 405 nm light).

Also chemically treated objects like fluorescing powder for the lifting of fingerprints or liquids for the indication of material damage (liquid penetrant testing) are visualized with blue light illumination and the use of orange filter goggles.

When using the 365 and 390 nm LED-Head, it is obligatory to wear transparent safety goggles (UV-shield). The 405 and 445 nm LED-Heads should be used with yellow or orange filter goggles as safety measure.

---

## 6. ACCESSORIES AND SPARE PARTS

---

9200.0002	Filter Goggles, Yellow (Long Pass Filter)
9200.0003	Filter Goggles, Orange (Long Pass Filter)
9200.0004	Filter Goggles, Dark Red (Long Pass Filter)
9200.0005	Filter Goggles*, Transparent (UV Protection)
9200.0050	Filter Goggles, Light Red (Long Pass Filter)
9200.0109	Filter Goggles / Green (Infrared and UV-protection)
9200.0007	Camera Filter, Dark Yellow
9200.0008	Camera Filter, Medium Yellow
9200.0009	Camera Filter, Orange
9200.0010	Camera Filter, Red
9200.0011	Camera Filter, Light Red
9200.0030	Case for Superlite M Sets
XXXX.XXXX	Liquid Light Guide Ø5 mm

---

\*Compliant with EN 166 and EN 170

## LED-Heads and Spare Parts

2350.0061	<b>M 05 / Charger</b>
2350.0016	<b>M 05 / Charger Adapter (only needed for Power Pack I)</b>
2350.0017	<b>M 05 / Power Pack I</b>
2350.0060	<b>M 05 / Power Pack II</b>
2350.0003	<b>M 05 / Lense UV (365 nm)</b>
2350.0002	<b>M 05 / Lense UV (390 nm)</b>
2350.0005	<b>M 05 / Lense Violet (405 nm)</b>
2350.0024	<b>M 05 / Lense Blue (445 nm)</b>
2350.0065	<b>M 05 / Lense Green (500 nm)</b>
2350.0001	<b>M 05 / Lense Green (530 nm)</b>
2350.0030	<b>M 05 / Lense yellow (595 nm)</b>
2350.0031	<b>M 05 / Lense red (620 nm)</b>
2350.0025	<b>M 05 / Lense IR (850 nm)</b>
2350.0062	<b>M 05 / Lense IR (940 nm)</b>
2350.0006	<b>M 05 / Lense VIS (400-700 nm)</b>
9200.0113	<b>M 05 / Wide angle attachment (only for Single-Wavelength-Heads)</b>
2350.0007	<b>M 05 / Lense Ultrathrow (6,500 K)</b>
2350.0026	<b>M 05 / LLG Adapter (LuD) VIS (6,500 K)</b>
2350.0027	<b>M 05 / LLG Adapter (LuD) Blue (445 nm)</b>
2350.0028	<b>M 05 / LLG Adapter (LuD) Green (530 nm)</b>
2350.0029	<b>M 05 / LLG Adapter (LuD) UV (365 nm)</b>
2350.0051	<b>M 05 / Multihead x6 (365 nm, 390 nm, 405 nm, 445 nm, 530 nm, 400-700 nm)</b>
2350.0052	<b>M 05 / Multihead x5 (365 nm, 390 nm, 405 nm, 445 nm, 400-700 nm)</b>
2350.0053	<b>M 05 / Multihead x4 (365 nm, 405 nm, 445 nm, 620 nm)</b>

## 7. INFORMATION ON UV RADIATION

The ultraviolet spectrum is divided into three regions:

UVA: 315-400 nm

UVB: 280-315 nm

UVC: 100-280 nm

UVA is a part of natural sunlight and its most visible effect on human tissue is the tan. As it is not visible, its intensity can be easily underestimated. An overdose of UVA will not cause UV-erythema (sunburn) but heat-erythema, like visible light. UVA can also cause biological long-term damages to human tissue.

For safety reasons it must absolutely be avoided to expose the unprotected eye or skin to ultraviolet radiation (see Safety Instructions, section 2).

UVB and UVC are more aggressive than UVA and can damage human tissue. However, the SUPERLITE **M 05** does not emit this radiation.

## 8. SPECTRAL RANGES AND APPLICATIONS

The following table shows spectral ranges and corresponding prominent applications in forensics.

Please note, that too intensive radiation may destroy a specimen. Start all examinations at low intensities.

Wavelength [nm]	Color	Goggles	Application
400 - 700	White	none	General search, footprints
320 - 500	UVA - Blue	Orange	General search
400 - 500	Blue	Orange	General search, sperm, urine, fibers etc.
320 - 400	UVA	Transparent (UV protection)	General search, stains, fingerprints, textile fibers, body fluids
400 - 420	Violet	none	Blood, finger marks in blood
		Yellow / Orange	Sperm, gun residuals
420 - 450	Indigo	Yellow / Orange	Sperm, finger prints treated with fluorescing powder
450 - 490	Blue	Orange	Fingerprints (Cyanoacrylate - Ardrox or Cyanoacrylate DFO)
490 - 520	Cyan	Orange	Fingerprints (Cyanoacrylate - Ninhydrin or Cyanoacrylate Rhodamine 6G)
520 - 560	Green	Orange	Fingerprints treated with DFO, reduction of background fluorescence
560 - 590	Yellow	Red	Fingerprints treated with Ninhydrin, reduction of background fluorescence
850 - 940	Infrared	IR-Camera	Gunshot residue and blood on dark background

All data are general guidelines only

---

## 9. TECHNICAL DATA

Model:	SUPERLITE M 05
LED type:	15 - 30 W
LED service life:	up to 15,000 hrs
Dimensions (L / Ø):	240 mm / 50 mm
Weight:	500 g
Power Pack:	7.2 V, 3 Ah Lithium-Ion
Operating temperature:	+10 - +35 °C
Rel. humidity:	<75% (<30 °C) <50% (>30 °C)
Operating altitude:	max. 2,000 m
Storage temperature:	-10 - +60 °C

Model:	SUPERLITE Charge Adapter
Power supply (AC):	110-230 VAC 50-60 Hz
Power supply (DC):	12 VDC
Power consumption:	20 W
Output:	5 V, 4 A
Operating temperature:	+10 - +35 °C
Rel. humidity:	<75% (<30 °C) <50% (>30 °C)
Operating altitude:	max. 2,000 m
Storage temperature:	-10 - +60 °C

Compliant with EN 55015, EN 61547 and  
EN 62471 Class 1

### Declaration of conformity

This unit conforms to all applicable EC Directives and corresponding harmonised standards. A written declaration of conformity can be supplied on request.

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LUMATEC GmbH