LUMATEC

SUPERLITE
S 04

THE PORTABLE LIGHT
SOURCE FOR
DIVERSE SPECTRA

Operating Instructions

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[EN] These operating instructions belong with this product. They contain important information for putting it into service and operating it. This should be noted also when this product is passed on to a third party.

Therefore look after these operating instructions for future reference!

A list of contents with the corresponding paged numbers can be found in the index on page 1.

LUMATEC prepared these instructions for experienced process engineers, technicians and manufacturing personnel. If you are new to UV light applications, do not understand the instructions or have any further questions and concerns, you can find our contact information on the last page of this booklet.

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Part No

2254.0000	SUPERLITE S 04	[Lamp: 120W AC; Filter Wheel: 10-fold, standard; Power Pack: Included]
2255.0000	SUPERLITE S 04	[Lamp: 120W AC; Filter Wheel: 10-fold, standard; Power Pack: Upgradeable]
2XXX.XXXX	SUPERLITE S 04	[Lamp / Filter Wheel / Battery: Customized version]

TABLE OF CONTENTS

1.	Intended Use.....	2
2.	Safety Instructions.....	2
3.	Equipment Overview	3
4.	Equipment Installation	4
5.	Equipment Operation	4
5.1	Light Source	4
5.2	Power Pack.....	5
6.	Special Instructions	7
7.	Equipment Maintenance and Service	8
8.	Accessories and Spare Parts	9
9.	Information on UV Radiation	9
10.	Spectral Ranges and Applications.....	10
11.	Technical Data.....	11

1. INTENDED USE

The light source SUPERLITE S 04, is solely intended for Non-Destructive Testing, fluorescence excitation and curing of polymers in industrial environments, such as work shops, (automated) assembly lines and laboratories.

2. SAFETY INSTRUCTIONS

Before using the SUPERLITE S 04, please read the following chapters of safety, installation and operation instructions.

Safety Symbols

Warning!

Potential Danger

Warning!

Optical Radiation

Mandatory Action!

Wear UV protection goggles

Mains Supply

The intelligent electronic power supply of the unit will adapt automatically to all AC currents from 110 to 240 Volts and 50 to 60 Hz. Connect power cord to a properly grounded AC outlet and make sure it is easily reachable during operation.

Explosive Surroundings

This unit is not meant for operation in explosive surroundings.

Humid Surroundings

This unit is only intended for use in dry environments.

Protection against Strong Optical Radiation

Never look directly at the light-emitting end of the light guide. Shield the operator or the light source. To avoid tissue damage, do not expose the unprotected eye or skin to the ultraviolet light. When work under unshielded UV radiation is necessary, wear UV protection goggles, long-sleeved shirt or lab coat and gloves.

Explosion of Mercury Vapor Lamps

In very rare occasions the lamp may burst and set free mercury to the environment. Remove all personnel from the room and ventilate thoroughly for 30 minutes. Any remains of mercury in the unit must be removed with mercury absorbent agent.

Danger of Fire

Do not deposit the light guide or connected optical equipment on inflammable objects like cloth or paper. The intensive radiation may be sufficient to set these objects on fire.

3. EQUIPMENT OVERVIEW


The SUPERLITE S 04 is a compact universal light source for the collection of evidence in the laboratory and at the scene of crime. By means of a filter wheel the user selects the spectrum appropriate for the individual job.

Light is delivered via a flexible light guide and released via a switch only when needed.

An optional battery as well as a carrying case enable independent operation at the scene of crime in dry locations.


Front view

- Light guide clip
- Stop screw
- Ball joint
- Telescopic arm
- Handle
- Shutter / Status light
- Intensity control knob
- Filter selection knob
- Light guide socket
- Power pack retaining clip
- Power pack level button
- Power pack status indicator


Rear view

- Name plate
- Key switch
- Charger socket
- Ventilation exhaust
- Mains fuses
- Mains socket
- Power pack holding clamp
- Power pack

- Collimator optics
- Light guide handle
- Sliding switch
- Safety clamp

4. EQUIPMENT INSTALLATION

Ventilation

Cooling air is let in at the bottom, while hot air is expelled through an opening at the rear of the unit. Both openings may not be covered in order to allow free ventilation. Minimum clearance should be 10-15 cm.

Obstruction of these openings will cause the unit to overheat.

UV Light Guide


Fully insert the flexible LIQUID LIGHT GUIDE with the Ø16 mm fitting into the light guide socket of the front panel of the light source. Treat the light guide with care, as it is an optical instrument. Do not kink or crush it and do not bend it too sharply to avoid light losses.

5. EQUIPMENT OPERATION

5.1 Light Source

Switching on

Turn on the power supply through the key switch on the back panel.

Off/Charge: Lamp is off, battery may be charged. Key can be removed. Posi-

tion for storage and transport.

Line/Charge: Lamp is on in mains mode, battery may be charged simultaneously. Key can not be removed.

Batt.: Lamp is on in battery mode. Key can not be removed.

After approximately three minutes the lamp reaches its rated power and is ready to operate.

Wavelength Selection

With the filter selection knob on the front panel you can choose between 10 different spectra. Emitted radiation is very high at the positions 400-700 nm, 300-500 nm and 400-500 nm. For safety reasons you must pull out the filter selection knob before selecting or deselecting these spectra. At short distances the illuminated area can become extremely bright, which may be harmful to the eye at extensive use. Turn down the intensity or make sure to take regular observation breaks.

Intensity Control

Use the intensity control knob to adjust the light intensity continuously, without changing the color temperature.

Releasing Light

Pressing the "SHUTTER" button on the front panel opens the shutter and the button will light up. Pressing it again closes the shutter.

Using the optional magnetic hand switch (LuD light guide socket only):

Switch locked

Light off

Light on

The shutter may also be opened by pulling and holding back the magnetic hand switch. Releasing the hand switch closes the shutter again. To avoid unintended exposure of light, the hand switch should be in locked position.

Telescopic Arm with Ball Joint and Clip

The telescopic arm may be used either as a stand to fix the light guide at an elevated position, or as an extended arm to guide the light guide, e.g. along the floor. To use it as a stand, mount the arm into the housing top and snap the light guide handle into the clip. To use it as an extended arm, snap the clip onto the light guide handle. Pull the telescopic arm to desired length and collapse it with small jolts. Use the ball joint to adjust the light guide to a desired position.

Status Light

The status light on the front panel of the unit can turn either white or red. White indicates, that the lamp is on. Red indicates, that the lamp has not ignited while switching on the unit, or may have shut down due to low battery level. Switch off the light source, as the status

indicator turns red. If the lamp does not ignite, exceeded lamp service life is mostly the cause. Lamp may burst! Replace lamp!

Switching off

The key switch will power off the lamp and separate the unit from the mains. Turn the key to the vertical position and remove it to prevent unintended or unauthorized operation of the light source.

After switching off the unit, it is recommended to wait three minutes before reigniting the lamp.

5.2 Power Pack

Charging

A fully charged SUPERLITE Power Pack will allow for approximately 90 minutes of independent operation. Charging time is four hours.

To check the charge condition, press the "LEVEL" button on the front panel of the power pack. Red segments on the status display indicate that the power pack should be recharged.

The power pack can be charged directly, or through the connected light source. Plug the provided charging cable into the corresponding socket on top of the power pack.

Use the “CHARGER” socket on the back panel of the light source, if you want to charge an attached power pack without unmounting it.

For charging, use the SUPERLITE Charger only! There are two options regarding power supply: 110-240 VAC and alternatively 12-24 VDC. Use the provided power cord to connect the charger with the mains supply. To connect to a 12-24 VDC power supply, use the “12-24 VDC” socket of the charger and the provided 12-24 VDC cable.


A green status LED of the charger indicates, that charging is in progress. After charging is completed the LED turns off.



If the power pack is attached to the light source during charging, the light source can be either off, or operate in mains mode. The battery can not be charged, if the light source is in battery mode. The light source can not be operated via the charger.

Mounting / Unmounting

The power pack has a rigid holding clamp in the rear and a spring loaded retaining clip with electrical connectors in the front.

Mount and unmount the power pack to the light source as per following instructions:

6. SPECIAL INSTRUCTIONS

Lamp Service Life

Assuming an average duty cycle of eight hours, a minimum lifetime of 2,500 hours can be expected for the installed lamp type. The lamp should not be ignited more often than necessary, as each ignition shortens lamp life. We recommend not to turn off the lamp, if it is needed again within the next three hours. As the lamp has a high internal pressure, it is possible that the bulb bursts in very rare occasions, especially if the specified lamp service life is exceeded. In this case mercury is emitted (see Safety Instructions, section 2).

Water Condensation

Water condensation may cause malfunction, if the unit is brought from a cold to a warm environment. Do not switch on the light source, until it has adapted to the higher temperatures. Alternatively, you may switch the unit on and wait for it to warm up, before moving the light source to warmer surroundings.

Overheating

If the maximum operating temperature is exceeded, the lamp is automatically switched off (internal thermostatic

switch). After the unit has cooled down sufficiently, the lamp can be ignited again. Reasons for overheating may be obstruction of ventilation openings (do not block the air intake and exhaust openings of the unit), or ambient temperatures exceeding specified maximum values.

Vapors of Solvents

Vapors of fluorinated or chlorinated hydrocarbon solvents will corrode the lamp and filters, even in small concentrations. Avoid these solvents in the vicinity of the unit. If this is not an option, place the unit as high as possible, as vapours are heavier than air.

7. EQUIPMENT MAINTENANCE AND SERVICE

 Caution, high voltages are present inside unit. Always disconnect power cord and unmount power pack before opening the unit. All maintenance and service activities should be carried out by trained specialists or be referred to an authorised service facility.

Exchanging the Lamp Module

 The lamp module should only be exchanged after the light source has cooled down completely.

Use the provided Torx key to open the four screws of the top cover of the unit.

Before removing the top cover, careful-

ly lift up its front side and unplug the grounding wire from the inner back end.


1. Unplug the lamp connector.
2. Open the safety screw above the lamp.
3. Pull the lamp module carefully up and lift it out of the two retaining pins.

Install the new lamp module in vice versa order. Before reassembling the cover, make sure the grounding wire is properly reconnected.

Never touch the filters!

Do not loosen any screws of the lamp module!

Used lamps contain mercury and must be disposed as hazardous waste.

8. ACCESSORIES AND SPARE PARTS

Accessories

9200.0013	UV Quartz Collimating Optics 15°
9200.0006	UV Collimating Optics 40°
9200.0014	Line Illuminator 50 mm
9200.0032	LitePad (Transmitted Light Inspections)
9200.0021	Spectral Magnifier / Camera Filter Set Ø77 mm Yellow, Orange, Red / Camera Adapter
9200.0005	Filter Goggles, Transparent (UV Protection)
9200.0002	Filter Goggles, Yellow (Long Pass Filter)
9200.0003	Filter Goggles, Orange (Long Pass Filter)
9200.0004	Filter Goggles, Red (Long Pass Filter)
9200.0012	Lorgnette, Yellow to Red (Long Pass Filters)
9200.0016	Telescopic Arm
9200.0059	SUPERLITE Power Pack
9200.0060	SUPERLITE Charger
9200.0061	Power Cord CN
9200.0062	Power Cord EU
9200.0063	Power Cord GB
9200.0064	Power Cord US
9200.0001	Soft Shell Field Case
9200.0017	Hard Shell Field Case
4440.700	Liquid Light Guide Ø3 mm, 700mm, Gooseneck
4502.1500	Liquid Light Guide Ø3 mm, 1,500 mm
4581.1800	Liquid Light Guide Ø5 mm, 1,800 mm, Magnetic Hand Switch

Spare Parts

2900.0003 Lamp Module 120 W AC

9. INFORMATION ON UV RADIATION

The ultraviolet spectrum is divided into three regions:

UVA: 315-400 nm

UVB: 280-315 nm

UVC: 100-280 nm

UVA is a part of natural sunlight and its most visible effect on human tissue is the tan. As it is not visible, its intensity can be easily underestimated. An overdose of UVA will not cause UV-erythema (sunburn) but heat-erythema, like visible light. UVA can also cause biological long-term damages to human tissue.

For safety reasons it must absolutely be avoided to expose the unprotected eye or skin to ultraviolet radiation (see Safety Instructions, section 2).

On customers request every unit can be fitted or retrofitted with filters for UVB and UVC.

UVB is more aggressive than UVA and also a part of natural sunlight. UVB has the same negative effects on human tissue as natural sunlight, the most common being the UV-erythema (sunburn).

UVC has an even greater detrimental effect on human tissue than UVB.

10. SPECTRAL RANGES AND APPLICATIONS

The following table shows spectral ranges of the standard 10-fold filter wheel

and corresponding prominent applications.

Please note, that too intensive radiation may destroy a specimen. Start all examinations at low intensities.

Wavelength [nm]	Color	Power [W]	Goggles	Application
400-700	White	10.0	none	General search, footprints
320-500	UVA - Blue	6.6	Orange	General search
400-500	Blue	5.5	Orange	General search, sperm, urine, fibers etc.
320-400	UVA	1.9	Transparent (UV protection)	General search, stains, fingerprints, textile fibers, body fluids
415	Violet	1.9	none	Blood, finger marks in blood
			Yellow / Orange	Sperm, gun residuals
440	Indigo Blue	2.1	Yellow / Orange	Sperm, finger prints treated with fluorescing powder
460	Royal Blue	2.0	Orange	Fingerprints (Cyanoacrylate - Ardrex or Cyanoacrylate DFO)
490	Turquoise	1.3	Orange	Fingerprints (Cyanoacrylate - Ninhydrin or Cyanoacrylate Rhodamine 6G)
550	Green - Yellow	1.4	Orange	Fingerprints treated with DFO, reduction of background fluorescence
570	Yellow	2.2	Red	Fingerprints treated with Ninhydrin, reduction of background fluorescence

All data are general guidelines only

11. TECHNICAL DATA

Model:	SUPERLITE S 04
Power supply:	110-240 VAC 50-60 Hz
Power consumption:	150 W
Protection:	2 fuses 2.5A slow
Lamp type:	Mercury vapour
Lamp service life:	up to 2,500 hrs
Dimensions (WxHxD):	175 x 140 x 260 mm
Weight:	3.8 kg
Operating temperature:	+10 - +35 °C
Rel. humidity:	<75% (<30 °C) <50% (>30 °C)
Operating altitude:	max. 2,000 m
Storage temperature:	-10 - +60 °C

Model:	SUPERLITE Power Pack
Output:	34 VDC
Dimensions (WxHxD):	175 x 50 x 260 mm
Weight:	2.4 kg
Operating temperature:	+10 - +35 °C
Rel. humidity:	<75% (<30 °C) <50% (>30 °C)
Operating altitude:	max. 2,000 m
Storage temperature:	-10 - +60 °C

Model:	SUPERLITE Charger
Power supply (AC):	110-240 VAC 50-60 Hz
Power supply (DC):	12-24 VDC
Power consumption:	110 W
Output:	34 VDC, 2 A
Dimensions (WxHxD):	54 x 64 x 145 mm
Weight:	0.6 kg
Operating temperature:	+10 - +35 °C

Rel. humidity:	<75% (<30 °C) <50% (>30 °C)
Operating altitude:	max. 2,000 m
Storage temperature:	-10 - +60 °C
Compliant with EN 60598, EN 60335, EN 61326	

Declaration of conformity

This unit conforms to all applicable EC Directives and corresponding harmonised standards. A written declaration of conformity can be supplied on request.

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LUMATEC GmbH